

Product Data Sheet

SERION coupling reagents

Product description	SERION coupling reagent N07 Ready-to-use activation reagent for the modification of negatively charged surfaces (e.g. carboxylate microparticles).
Order No.	CRN07
Lot No.	Lot specific
Matrix	Buffer containing 0.02 %Vol. Tween™ 20 and 0.1 %Vol. ProClin™ 950
Preservatives	0.1 %Vol. ProClin™ 950
Appearance	Clear colorless solution
Properties	Shifts isoelectric point of carboxylate microparticles to lower pH values. Determined by zetapotential-pH-titration with MagSERION carboxy beads.
Storage	2 – 8 °C (shipment at room temperature) Do not freeze.
Date of manufacture	Lot specific
Expiry date	No expiry date has been assigned for this product. Users should determine stability in their own system.
Manufacturer	Institut Virion\Serion GmbH, Germany

For research or manufacturing use only. Not for use in diagnostic procedures.

Date: 2019-12-12