

## Product Data Sheet

### SERION coupling reagents

Product description	SERION coupling reagent P01 Ready-to-use activation reagent for the modification of positively charged surfaces (e.g. amine microparticles).
Order No.	CRP01
Lot No.	Lot specific
Matrix	Buffer containing 0.02 %Vol. Tween™ 20 and 0.1 %Vol. ProClin™ 950
Preservatives	0.1 %Vol. ProClin™ 950
Appearance	Clear colorless solution
Properties	Shifts isoelectric point of amine microparticles to lower pH values. Determined by zetapotential-pH-titration with MagSERION amine beads.
Storage	2 – 8 °C (shipment at room temperature) Do not freeze.
Date of manufacture	Lot specific
Expiry date	No expiry date has been assigned for this product. Users should determine stability in their own system.
Manufacturer	Institut Virion\Serion GmbH, Germany

*For research or manufacturing use only. Not for use in diagnostic procedures.*

Date: 2019-12-12