

Product Data Sheet

SERION blocking buffer for coupling reagents

Product description	SERION blocking buffer for coupling reagents Ready-to-use buffer for blocking of microparticles (e.g. MagSERION carboxy or amine beads) after biomolecule coupling.
Order No.	CB01
Lot No.	Lot specific
Preservatives	0.09 % sodium azide
pH-value	pH 8.4
Appearance	Clear colorless solution
Storage	2 – 8 °C (shipment at room temperature) Do not freeze.
Expiry date	36 months from date of manufacture (under recommended storage conditions)
Manufacturer	Institut Virion\Serion GmbH, Germany

For research or manufacturing use only. Not for use in diagnostic procedures.

Date: 2019-12-16